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## DABTunel assay report (Paraffin section)

### 1. Experimental principle

When cells undergo apoptosis, the double strands of chromosomal DNA break, generating a large number of sticky 3'-OH termini. Under the action of terminal deoxynucleotidyl transferase (TdT), derivatives formed by deoxynucleotides and fluorescein, peroxidase, alkaline phosphatase, or biotin can be labeled to the 3'-terminus of DNA, thereby enabling the detection of apoptotic cells (TUNEL method).

### 2. Laboratory equipment and reagents

#### 2.1 Laboratory equipment

Equipment	Manufacturers	Model
Dehydrator	Wuhan Junjie Electronics Co., Ltd	JT-12S
Paraffin embedding machine	Wuhan Junjie Electronics Co., Ltd	JB-P5
Frozen table	Wuhan Junjie Electronics Co., Ltd	JB-L5
Rotary microtome	Shanghai Leica Instrument Co., Ltd	HistoCoreBIOCUT
Tissue machine	Wuhan Junjie Electronics Co., Ltd	JK-5
Oven	Tianjin Leibo Terry Equipment Co., Ltd	GFL125
Glass microscope slides	Haimen Shenying Experimental Equipment Factory	188109
Rocker	Servicebio	SYC-Z100
Vortex	Servicebio	MX-F
Micro-centrifuge	Servicebio	DS-S 100

Microwave oven	Midea	M1-L213B
Liquid Blocker PAP Pen	Gene tech	GT1001
Pipettor	Dragon	KE003068
Coverslips	Jiangsu Huida Medical Instruments Co., Ltd	710510
Microscope	NIKON	ECLIPSE E100
Imaging System	Konfoong Bioinformation Tech Co.,Ltd	KF-PRO-120

## 2.2 Laboratory reagents

Reagents	Manufacturers	Catlog
Anhydrous ethanol	Hangzhou Hongda Chemical Instrument Co., LTD	SJ003614
Xylene	Sinopharm Chemical Reagent co., Ltd.	10023418
PBS solution	Hangzhou Haoke Biotechnology Co., Ltd	HK0002
Citric acid antigen repair solution (9.0)	Hangzhou Haoke Biotechnology Co., Ltd	HKI0004
DABtunel Assay Kit	Hangzhou Haoke Biotechnology Co., Ltd	HKI0012
Hematoxylin staining solution	Hangzhou Haoke Biotechnology Co., Ltd	HK2053
Hematoxylin differentiate solution	Hangzhou Haoke Biotechnology Co., Ltd	HK2054
Hematoxylin bluing solution	Hangzhou Haoke Biotechnology Co., Ltd	HK2055
Neutral gum	Sinopharm Chemical Reagent co., Ltd.	10004160

### 3. Experimental procedure

**3.1 Dewaxing and hydration:** Put the sections into xylene I 12 min - xylene II 12 min - anhydrous ethanol I 6 min - 95% Ethyl alcohol for 6 min - 85% Ethyl alcohol for 6 min, rinsing with tap water for 2 min.

**3.2 Repair:** The tissue sections were placed in a repair box filled with citric acid antigen repair solution (9.0), and the antigen repair was carried out in a microwave oven at medium heat for 8min. During this process, excessive evaporation of buffer solution should be prevented and the slices should not be dried. After natural cooling, the slide was placed in double steaming water and washed by shaking on the decolorizing table for 3 times, 5 min each time.

**3.3 Draw a circle around the tissue:** After the section is slightly dry, draw a small circle with a histochemical pen along the outer outline of the tissue with a distance of 3-4 mm from the tissue. After drawing, wash it with pure water for 3 times, 5 min each time.

**3.4 Add reagents:** Appropriate amount of reagent TdT enzyme and white light reaction solution were taken from tunel kit according to the number of slides and tissue size, mixed at 1:100, then added to the ring to cover the tissue, the slices were placed flat in a wet box, incubated at 37°C for 1 hour, and a small amount of water was added to the wet box to maintain humidity.

**3.5 HRP secondary antibody incubation:** Sections were washed 3 times with PBS (PH7.4) for 5min each time. Then, HRP secondary antibody working liquid was added and incubated at room temperature for 30 min, the working liquid concentration HRP: PBST=1:500, and the color was developed by DAB after washing.

**3.6 Restaining of nucleus:** Hematoxylin restained the nuclei for 2-3min. Wash with tap water, hematoxylin differentiate solution differentiation for 2s, rinse with tap water, hematoxylin bluing solution back blue 15-30s, water rinse.

**3.7 Dehydration and sealing:** Put the sections into 75% Ethyl alcohol for 4 min - 85% Ethyl alcohol for 4 min - Anhydrous ethanol I for 4 min - Anhydrous ethanol II for 4 min - xylene I for 4 min - xylene II for 4 min, after slightly drying, , take the slices out of xylene for a little dry, and seal them with neutral gum.

#### **4. Interpretation of results**

Hematoxylin stained nuclei are blue, and DAB positive apoptotic nuclei are brownish yellow

#### **5. Precautions**

- 5.1 Pay attention to whether the slice dewaxing is thorough;
- 5.2 Don't let the slices dry out during the experiment;
- 5.3 Wash with pure water before adding tunel reaction solution.