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Blood Routine Test Report

Experimental principle

Blood routine is the most general, the most basic blood test. Blood consists of two parts: liquid and visible cells. Blood routine tests the cellular part of blood. Blood has three different types of functional cells - red blood cells (commonly known as red blood cells), white blood cells (divided into five categories including lymphocytes, monocytes, neutrophils, eosinophils, basophilic granulocytes), and platelets. The disease was judged by observing the quantity change and morphological distribution. It is one of the commonly used auxiliary examination means for doctors to diagnose the condition. Semiconductor laser scattering technology combined with flow cytometry technology to achieve accurate five classification of white blood cells, environmentally friendly cyanide-free determination of hemoglobin.

Laboratory equipment and reagents

Laboratory equipment

Equipment	Manufacturers	Model
Five-part fully automated hematology analyzer	Tecom	TEK8500H4-0502

Laboratory reagents

Reagents	Manufacturers	Catlog
Sheath fluid for blood cell analysis	Tecom	A
Hemolytic agents for blood cell analysis	Tecom	TEK8510
Cleaning solution for blood cell analysis	Tecom	C1

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Operation steps

Boot → select animal type → select detection mode → sample loading detection

Result

See the results table for details

Sample and shipping requirements

For sample, take an appropriate amount of whole blood with an EDTA anticoagulant tube. The optimal concentration ratio of EDTA to blood is 1.5mg/mL. After sampling, gently turn up and down several times to mix the blood well with the anticoagulant.

For shipping conditions, the samples are shipped at 4°C. Please do not contact with ice packs, and do not violently impact.

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Biotechnology Co., Ltd.