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Real Time PCR Report

Experimental principle

Quantitative Real-time PCR is a method in which a fluorescent chemical is used to measure the total amount of product after each PCR cycle in DNA amplification. A method of quantitative analysis of a specific DNA sequence in a test sample by internal or external parameters.

Laboratory equipment and reagents

Laboratory equipment

Equipment	Manufacturers	Model
Tissue Grinder	Servicebio	KZ-II
High speed refrigerated centrifuge	Haier	LX-165T2R
Real time PCR System	ABI	7500
Clean bench	BIOBASE	BBS-DDC
Ultramicro spectrophotometer	Canhelp Genomics Co., Ltd	CANHELP-100
Laboratory pure water system	Zhiang Instrument (Shanghai) Co., Ltd	Clever-S

Laboratory reagents

Reagents	Manufacturers	Catlog
RNA Extraction	Hangzhou Haoke Biotechnology Co., Ltd	HKR022
trichloromethane	Sinopharm Chemical Reagent co., Ltd.	10006818

Isopropyl alcohol	Sinopharm Chemical Reagent co., Ltd.	80109218
Anhydrous ethanol	Sinopharm Chemical Reagent co., Ltd.	10009218
All-in-One First-Strand cDNA Synthesis SuperMix for qPCR	TRAN	AE341
2X SYBR Green qPCR Master Mix	APExBIO	K1070

Experimental method

Isolate RNA

- 1) Take homogenate tube, add 1mL of Trizol Reagent, place on ice for pre-cooling.
- 2) Take 100mg tissue and add it to the homogenate tube.
- 3) The homogenizer fully grinds until no tissue blocks are visible.
- 4) Centrifuge at 12000rpm for 10min to take the supernatant.
- 5) Add 250 μ L trichloromethane, reverse the centrifuge tube for 15s, mix thoroughly, and let stand for 3min.
- 6) Centrifuge at 13000rpm at 4 $^{\circ}$ C for 15min.
- 7) Transfer the supernatant to the purification tube, add 1.2 times the volume of anhydrous ethanol, and mix well.
- 8) Centrifuge at 12000rpm at 4 $^{\circ}$ C for 5min, and purify the white precipitate on the tube filter membrane as RNA.
- 9) Discard the liquid in the lower part of the purification tube, add 75% ethanol 650 μ L at 4 $^{\circ}$ C, centrifuge at 12000rpm for 5min and wash.
- 10) Repeat washing 2-3 times.
- 11) Centrifuge at 12000rpm at 4 $^{\circ}$ C for 5min, discard the residual liquid, place the centrifuge tube on a super-clean table and blow for 3-5min.
- 12) Add 50 μ L of RNA-free water to dissolve RNA.

13) Use an ultramicrospectrophotometer to detect RNA concentration and purity: After blank zero adjustment of the instrument, take 2 μ L RNA solution to be tested and put it on the detection base, lower the sample arm, and use the software on the computer to start the absorption value detection.

14) The RNA with excessive concentration is diluted in an appropriate proportion so that its final concentration is 100-500 ng/ μ L.

First Strand cDNA Synthesis

- 1) Take a PCR tube and add a solution containing 1 μ g RNA.
- 2) Add 1 μ L dsDNASE and 1 μ L 10*dsDNase Buffer.
- 3) Fill with ribonuclease free deionized water to 10 μ L and mix well.
- 4) Incubate at 37 $^{\circ}$ C for 2min on PCR apparatus to remove genomic DNA contamination. dsDNAase was inactivated after incubation at 65 $^{\circ}$ C for 2min.
- 5) Add 4 μ L All-in-One First-Strand Synthesis MasterMix.
- 6) Fill with ribonuclease free deionized water to 20 μ L and mix well.
- 7) Incubation at 50 $^{\circ}$ C for 15min on the PCR apparatus, and then incubation at 85 $^{\circ}$ C for 5min to inactivate reverse transcriptase.

Preparation of PCR Master Mix

- 1) Take 0.2mL PCR tubes and prepare the reaction system as follows, with 3 tubes for each retroproduct.

Reagent	Volume
cDNA	1 μ L
Forward Primer (10 μ M)	0.4 μ L
Reverse Primer (10 μ M)	0.4 μ L
Taq SYBR Green qPCR Premix	10 μ L

ddH₂O

To 20 μ L

2) PCR amplification

Stage1 predegeneration	Stage2 (50cycles)	Stage3 (Melt Curve)
94°C, 2min	94°C, 5s→58°C, 10s→72°C, 15s	60°C→95°C

Data analysis

$\Delta\Delta$ CT method:

A=CT(target gene, sample)- CT(internal standard gene, sample)

B=CT(target gene, control)- CT(internal standard gene, control)

K=A-B

RNA Expression = 2^{-K}