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Paraffin sections polychromatic immunofluorescence experiment report (TSA)

1. Experimental principle

Tyramide signal amplification (TSA) technique is a kind of enzymatic detection method using horseradish peroxidase (HRP) to label target protein, which is similar to the conventional immunohistochemical DAB color development method. TSA technology also uses HRP labeled secondary antibodies, and also has corresponding "color rendering" steps to catalyze the addition of tyramine fluorescein substrate in the reaction system to produce an activated fluorescent substrate, which can covalently bind tyramine fluorescein and other residues on the antigen to make the sample covalently bind tyramine fluorescein. After that, the non-covalently bound primary - secondary -HRP complex was washed by thermal repair method, and the next primary -hrp secondary - antibody was repeated for the second round of incubation, so that multiple labeling could be achieved repeatedly.

2. Laboratory equipment and reagents

2.1 Laboratory equipment

Equipment	Manufacturers	Model
Dehydrator	Wuhan Junjie Electronics Co., Ltd	JT-12S
Paraffin embedding machine	Wuhan Junjie Electronics Co., Ltd	JB-P5
Rotary microtome	Shanghai Leica Instrument Co., Ltd	HistoCoreBIOCUT
Frozen table	Wuhan Junjie Electronics Co., Ltd	JB-L5

Tissue machine	Wuhan Junjie Electronics Co., Ltd	JK-5
Oven	Tianjin Leibo Terry Equipment Co., Ltd	GFL125
Coverslips	Jiangsu Huida Medical Instruments Co., Ltd	710510
Glass microscope slides	Haimen Shenying Experimental Equipment Factory	188109
Microwave oven	Midea	M1-L213B
Rocker	Servicebio	SYC-Z100
Vortex	Servicebio	MX-F2
Micro-centrifuge	Servicebio	DS-S 100
Pipettor	Dragon	KE003068
Liquid Blocker PAP Pen	Gene tech	GT1001
Fluorescence microscope	NIKON	ECLIPSE C1
3D Imaging System	3D HISTECH	Pannoramic SCAN

2.2 Laboratory reagents

Reagents	Manufacturers	Catlog
Anhydrous ethanol	Hangzhou Hongda Chemical Instrument Co., Ltd	SJ003614
Xylene	Sinopharm Chemical Reagent Co., Ltd	10023418
EDTA Antigen repair solution (PH8.0)	Hangzhou Haoke Biotechnology Co., Ltd	HKI0003
PBS solution	Hangzhou Haoke Biotechnology Co., Ltd	HK0002

BSA	Hangzhou Haoke Biotechnology Co., Ltd	HKW2084
Primary antibody: xxxx		
Second antibody: HRP	Hangzhou Haoke Biotechnology Co., Ltd	
hypersensitive goat anti-rabbit secondary antibody		HKI0029
Flare520	Hangzhou Haoke Biotechnology Co., Ltd	HKI0014
DAPI	Hangzhou Haoke Biotechnology Co., Ltd	HKI0005
Anti-fluorescence quenching sealing tablets	Hangzhou Haoke Biotechnology Co., Ltd	HKI0007-1

3. Experimental procedure

3.1 Dewaxing and hydration: Put the sections into xylene I 12 min - xylene II 12 min - anhydrous ethanol I 6 min - 95% Ethyl alcohol for 6 min - 85% Ethyl alcohol for 6 min, rinsing with distilled water.

3.2 Repair: The tissue sections were placed in a repair box filled with EDTA repair solution (PH8.0) and the antigen repair was carried out in a microwave oven. The repair procedure was: medium fire for 8min, ceasefire for 8min, and medium low fire for 7min. During this process, excessive evaporation of buffer solution should be prevented and the pieces should not be dried. After natural cooling, the slide was placed in PBS (PH7.4) and washed by shaking on the decolorizing shaker for 3 times, 5min each time.

3.3 Block endogenous peroxidase : The sections were incubated with 3% hydrogen peroxide at room temperature for 25min. The slides were placed in PBS (PH7.4) and washed three times on a decolorizing shaker for 5min each time.

3.4 Draw a circle around the tissue: Use a histochemical pen to draw a small circle 3-4 mm apart from the tissue along the outer outline of the tissue, and then add enough PBS to ensure that the subsequent addition of blocking serum, primary antibody, secondary

antibody, and color development agent can completely cover the tissue without flowing along the sections.

3.5 Serum blocking: After the section is slightly dried, draw a circle around the tissue with a histochemical pen (to prevent the antibody from flowing away), and add 3%BSA to the circle and seal it at room temperature for more than 30 minutes.

3.6 Add primary antibody: After the section is slightly dried, draw a circle around the tissue with a tissue pen (to prevent the antibody from flowing away), and add a certain proportion of dilute antibody to cover the tissue in the circle. The slices were incubated flat in a wet box at 4°C overnight.

3.7 Add secondary antibody: The sections were placed in PBS (PH7.4) and washed by shaking on the decolorizing shaker for 3 times, 5min each time. After the slices were slightly dried, the tissue covered by the second antibody of the corresponding species of the first antibody was added to the ring and incubated at room temperature for 50min away from light.

3.8 TSA: The sections were placed in PBS (PH7.4) and washed by shaking on the decolorizing shaker for 3 times, 5min each time. Add the corresponding color Flare signal amplification reagent, 3-5min. Place the slide in PBS (PH7.4), shake and wash it on the decolorizing shaker for 3 times, 5min each time. After this step, repeat steps 2-7, label the second or third antibody, and then proceed to step 8.

3.9 DAPI restaining nuclei: The sections were placed in PBS (PH7.4) and washed by shaking on the decolorizing table for 3 times, 5 min each time. Blot the excess PBS in the ring and then add DAPI dye solution, and incubate at room temperature for 8min away from light.

3.10 Sealing: The sections were placed in PBS (PH7.4) and washed by shaking on the decolorizing table for 3 times, 5min each time. Remove the slices and seal them with anti-fluorescence quenching tablets.

3.11 Microscopy photography: The sections were observed and the images were collected under a fluorescence microscope, and the corresponding channels were selected according to the corresponding wavelength.

4. Interpretation of results

The nuclei stained by DAPI are blue under ultraviolet excitation, and the positive expression is the color labeled by the corresponding fluorescent reagent (refer to the table below for specific fluorescence channels).

Flare Dye	Excitation wavelength	Emission wavelength
DAPI	350	420
Flare480	450	480
Flare520	490	520
Flare570	550	570
Flare620	590	620
Flare690	630	690
Flare780	750	780

5. Precautions

- 5.1 Pay attention to whether the slice dewaxing is thorough;
- 5.2 Don't let the slices dry out during the experiment;
- 5.3 The addition of fluorescence amplification reagent is irreversible and should be added with caution.